

# ***Stand Perimeter Schweigger Type User's Manual***



3700 Sagamore Parkway North  
Lafayette, IN 47904 USA  
Tel: (765) 423-1505 • (800) 428-7545  
Fax: (765) 423-4111  
[info@lafayetteinstrument.com](mailto:info@lafayetteinstrument.com)  
[www.lafayetteinstrument.com](http://www.lafayetteinstrument.com)

## Description

For optometric, clinical and laboratory experimentation work. The Perimeter is used to map the areas of the retina sensitive to stimuli of different colors. Included with the apparatus are stimulus patches of red, yellow, green, blue, and white, 1/8" square, on a 5/8" aluminum disk.

The Perimeter is constructed in such a manner that measurements can be made of the threshold between regions where the stimulus is and is not discriminated. This is accomplished by having Subject (S) fixate on central spot (small mirror) while Experimenter (E) introduces the stimulus patch by moving it along the arm of the Perimeter. The distance of the stimulus from the eye is held constant by moving the stimulus target along the slit opening of the curved metal arm of the Perimeter, the eye being located at the center of the curvature of the arm. By moving the stimulus along the arm, E establishes the boundary between the regions of successful and unsuccessful color discrimination. The spatial position of the threshold is expressed by the angle between the line of fixation and the line from the eye to the threshold. Furthermore, by rotating the Perimeter arm about its pivot, threshold measurements along any line radiating out from the center of the visual field can be made; and by moving the stimulus patch inward and outward, threshold measurements for direction of stimulus movement can be made.

## Procedure

Give S a comfortable seat and place cheekbone against back plastic support so that his left eye is directly centered in the fixation mirror. Have S keep his right eye covered during the experiment.

The nine rows of Table 2 show the angular positions of the Perimeter arm at which measurements are to be made. The columns give the hues of the stimulus patches and specify the direction in which E must move the stimulus.

In order that S be as naive as possible with regard to the color being used, employ the following table or randomized stimulus sequences.

**Table 1: A Table for the Experimenter to Establish the Sequence of Colors**

<u>Row</u>	<u>Sequence</u>	<u>Row</u>	<u>Sequence</u>
1	GYBWR	5	GYRBW
2	BRWGY	6	RBGWY
3	YBRWG	7	BWRYG
4	YWGBR	8	RGYWB

Choose a row from Table 1—say Row 3—with the arm set at 0 degrees, present the stimuli in the order shown i.e., yellow, blue, red, white, and green. Present each color twice—first moving inward, then outward, beginning with the stimulus patch at the periphery and slowly moving the patch into the field of vision. The S must have no foreknowledge of the color for the in procedure. Tell S to report aloud and continuously what he sees—for example, “nothing, nothing, gray, blue, blue, yellow.” When the color is identified, stop and then read and record the angular distance from the center of fixation under in and yellow in the first row of Table 2. Next, administer the out trial, starting well within S’s field of color vision, and slowly moving the yellow patch toward the periphery. Again instruct S to make a continuous oral report. Do not stop movement until S reports a change in the color.

After completing the above procedure for the five colors, rotate the metal arm 45 degrees to obtain the data required for the second row of Table 2. The sequence of colors may again be determined by a random selection from Table 1 or you may simply proceed to the next available row returning to Row 1 following the completion of Row 8.

**Table 2: Data Sheet**

Angle of Rotation	Red		Yellow		Green		Blue		White	
	In	Out	In	Out	In	Out	In	Out	In	Out
0										
45										
90										
135										
180										
225										
270										
315										
360										

# Terms and Conditions

## LIC Worldwide Headquarters

Toll-Free: (800) 428-7545 (USA only)

Phone: (765) 423-1505

Fax: (765) 423-4111

Email: [sales@lafayetteinstrument.com](mailto:sales@lafayetteinstrument.com)

[export@lafayetteinstrument.com](mailto:export@lafayetteinstrument.com) (Outside the USA)

## Mailing Address:

Lafayette Instrument Company

PO Box 5729

Lafayette, IN 47903, USA

## Lafayette Instrument Europe:

Phone: +44 1509 817700

Fax: +44 1509 817701

Email: [esusales@lafayetteinstrument.com](mailto:esusales@lafayetteinstrument.com)

## Phone, Fax, Email or Mail-in Orders

All orders need to be accompanied by a hard copy of your purchase order. All orders must include the following information:

- Quantity
- Part Number
- Description
- Your purchase order number or method of pre-payment
- Your tax status (include tax-exempt numbers)
- Shipping address for this order
- Billing address for the invoice we'll mail when this order is shipped
- Signature and typed name of person authorized to order these products
- Your telephone number
- Your email address
- Your FAX number

## Domestic Terms

There is a \$50 minimum order. Open accounts can be extended to most recognized businesses. Net amount due 30 days from the date of shipment unless otherwise specified by us. Enclose payment with the order; charge with VISA, MasterCard, American Express, or pay COD. We must have a hard copy of your purchase order by mail, E-mail or fax. Students, individuals and private companies may call for a credit application.

## International Payment Information

There is a \$50 minimum order. Payment must be made in advance by: draft drawn on a major US bank; wire transfers to our account; charge with VISA, MasterCard, American Express, or confirmed irrevocable letter of credit. Proforma invoices will be provided upon request.

## Exports

If ordering instrumentation for use outside the USA, please specify the country of ultimate destination, as well as the power requirements (110V/60Hz or 220V/50Hz). Some model numbers for 220V/50Hz will have a "C" suffix.

## Quotations

Quotations are supplied upon request. Written quotations will include the price of goods, cost of shipping and handling, if requested, and estimated delivery time frame. Quotations are good for 30 days, unless otherwise noted. Following that time, prices are subject to change and will be re-quoted at your request.

## Cancellations

Orders for custom products, custom assemblies or instruments built to customer specifications will be subject to a cancellation penalty of 100%. Payment for up to 100% of the invoice value of custom products may be required in advance. Cancellation for a standard Lafayette Instrument manufactured product once the product has been shipped will normally be assessed a charge of 25% of the invoice value, plus shipping charges. Resell items, like custom products, will be subject to a cancellation penalty of 100%.

## Exchanges and Refunds

Please see the cancellation penalty as described above. No item may be returned without prior authorization of Lafayette Instrument Company and a Return Goods Authorization (RGA#) number which must be affixed to the shipping label of the returned goods. The merchandise should be packed well, insured for the full value and returned along with a cover letter explaining the reason for return. Unopened merchandise may be returned prepaid within thirty (30) days after receipt of the item and in the original shipping carton. Collect shipments will not be accepted. Product must be returned in saleable condition, and credit is subject to inspection of the merchandise.

## >Returns

Instrumentation may not be returned without first receiving a Return Goods Authorization Number (RGA). When returning instrumentation for service,

please call Lafayette Instrument to receive a RGA number. Your RGA number will be good for 30 days. Address the shipment to:

Lafayette Instrument Company

3700 Sagamore Parkway North

Lafayette, IN 47904, USA.

Shipments cannot be received at the PO Box. The items should be packed well, insured for full value, and returned along with a cover letter explaining the malfunction. An estimate of repair will be given prior to completion ONLY if requested in your enclosed cover letter. We must have a hard copy of your purchase order by mail or fax, or repair work cannot commence for non-warranty repairs.

## Damaged Goods

Damaged instrumentation should not be returned to Lafayette Instrument prior to a thorough inspection. If a shipment arrives damaged, note damage on delivery bill and have the driver sign it to acknowledge the damage. Contact the delivery service, and they will file an insurance claim. If damage is not detected at the time of delivery, contact the carrier/shipper and request an inspection within 10 days of the original delivery. Please call the Lafayette Instrument Customer Service Department for repair or replacement of the damaged merchandise.

## Limited Warranty

Lafayette Instrument Company warrants equipment manufactured by the company to be free of defects in material and workmanship for a period of one year from the date of shipment, except as provided hereinafter. The original manufacturer's warranty will be honored by Lafayette Instrument for items not manufactured by Lafayette Instrument Company, i.e. resell items. This assumes normal usage under commonly accepted operating parameters and excludes consumable products.

Warranty period for repairs or used instrumentation purchased from Lafayette Instrument is 90 days. Lafayette Instrument Company agrees either to repair or replace, at its sole option and free of part charges to the customer, instrumentation which, under proper and normal conditions of use, proves to be defective within the warranty period. Warranty for any parts of such repaired or replaced instrumentation shall be covered under the same limited warranty and shall have a warranty period of 90 days from the date of shipment or the remainder of the original warranty period whichever is greater. This warranty and remedy are given expressly and in lieu of all other warranties, expressed or implied, of merchantability or fitness for a particular purpose and constitutes the only warranty made by Lafayette Instrument Company.

Lafayette Instrument Company neither assumes nor authorizes any person to assume for it any other liability in connection with the sale, installation, service or use of its instrumentation. Lafayette Instrument Company shall have no liability whatsoever for special, consequential, or punitive damages of any kind from any cause arising out of the sale, installation, service or use of its instrumentation. All products manufactured by Lafayette Instrument Company are tested and inspected prior to shipment. Upon prompt notification by the Customer, Lafayette Instrument Company will correct any defect in warranted equipment of its manufacture either, at its option, by return of the item to the factory, or shipment of a repaired or replacement part. Lafayette Instrument Company will not be obliged, however, to replace or repair any piece of equipment, which has been abused, improperly installed, altered, damaged, or repaired by others. Defects in equipment do not include decomposition, wear, or damage by chemical action or corrosion, or damage incurred during shipment.

## Limited Obligations Covered by this Warranty

1. In the case of instruments not of Lafayette Instrument Company manufacture, the original manufacturer's warranty applies.
2. Shipping charges under warranty are covered only in one direction. The customer is responsible for shipping charges to the factory if return of the part is required.
3. This warranty does not cover damage to components due to improper installation by the customer.
4. Consumable and/or expendable items, including but not limited to electrodes, lights, batteries, fuses, O-rings, gaskets, and tubing, are excluded from warranty.
5. Failure by the customer to perform normal and reasonable maintenance on instruments will void warranty claims.
6. If the original invoice for the instrument is issued to a company that is not the company of the end user, and not an authorized Lafayette Instrument Company distributor, then all requests for warranty must be processed through the company that sold the product to the end user, and not directly to Lafayette Instrument Company.

## Export License

The U.S. Department of Commerce requires an export license for any polygraph system shipment with an ULTIMATE destination other than: Australia, Japan, New Zealand or any NATO Member Countries. It is against U.S. law to ship a Polygraph system to any other country without an export license. If the ultimate destination is not one of the above listed countries, contact us for the required license application forms.