

Polygraph Pro Suite v 4.6

v 4.6.3.2

Released: November 5, 2024

Patches

- Resolved issue that file created by “Export Scientific Measures” did not contain useful information
- Resolved issue that Menu item for Updating software opened link that is no longer valid
- Modified NCCA export to conform to standards
 - Swapped data contained in EDA1 and EDA2 columns to match LXEdge
 - EDA1 now contains “EDA RAW” sensor data and will be imported into LXEdge as “Manual”
 - EDA2 now contains “EDA” sensor data and will be imported into LXEdge as “Automatic”
 - Resolved filenames being incorrect
 - Events list is now properly implemented
 - Fixed certain cases where valid export data could be negative
 - Added missing start and stop questions ("X" and "XX")
 - Added label handling for DLST exams; these exams no longer repeat the question type
 - Modified question set templates to adhere to NCCA output standards
 - Removed ability to edit or add question types
 - Show the correct value “Number of channels”
 - Resolved issue with NCCA exporting masked PFs

v 4.6.2.1

Released: August 21, 2024



Polygraph Pro Suite Release Notes

Patches

- Resolved issue that Menu item for Updating software opened link that is no longer valid
- Modified NCCA export to conform to standards
 - Swapped data contained in EDA1 and EDA2 columns to match LXEdge
 - EDA1 now contains “EDA RAW” sensor data and will be imported into LXEdge as “Manual”
 - EDA2 now contains “EDA” sensor data and will be imported into LXEdge as “Automatic”
 - Resolved filenames being incorrect
 - Events list is now properly implemented
 - Fixed certain cases where valid export data could be negative
 - Added missing start and stop questions ("X" and "XX")
 - Added label handling for DLST exams; these exams no longer repeat the question type
 - Modified question set templates to adhere to NCCA output standards
 - Removed ability to edit or add question types
 - Show the correct value “Number of channels”
 - Resolved issue with NCCA exporting masked PFs

v 4.6.1.1

Released: June 3, 2024

Patches

- Resolved intermittent errors when loading the scoring application.

v 4.6.1

Released: April 26, 2024

Patches

- Resolved issues when loading settings so that mixed culture files will correctly be interpreted.
- Resolved issues when saving settings so that decimal values are always written in Invariant Culture format, thus preventing mixed culture files.

v 4.6.0.3

Released: November 29, 2023

Patches

- Updated installer to register necessary DLL with Windows.

v 4.6.0.2

Released: August 14, 2023

Patches

- Resolved parsing error of Chart Recorder when reading Calibrate.ini

Polygraph Pro Suite Release Notes

v 4.6.0.1

Released: August 11, 2023

Patches

- Implemented protections to cause Settings regeneration to be multi-threading safe.

v 4.6.0.0

Released: August 10, 2023

Patches

- Resolve issue where Chart Viewer would crash when trying to read a setting value that does not exist.